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<b>Code</b>	<b>Description and Possible Remedy</b>
2353	<p><b>Expected plug/hole/filter not found in filter cube</b></p> <p>This error indicates that the filter cube is not installed and is required for the read. Install the filter cube or check that it is installed correctly.</p>
2700	<p><b>Error attempting to run the barcode scanner SET command</b></p> <p>A response returned from the scanner is invalid.</p>
2701	<p><b>Error attempting to run the barcode scanner SET command</b></p> <p>The command message is calling out an invalid barcode location. Valid numbers are 1–4.</p>
2702	<p><b>Error attempting to get barcode scanner information regarding one of the barcode types</b></p> <p>The barcode type returned by the scanner is not one of those expected.</p>
2703	<p><b>Barcode type is not supported</b></p> <p>One of the four barcode types is not supported by the scanner.</p>
2704	<p><b>Error disabling start/stop character transmission</b></p> <p>While attempting to tell the scanner to disable the transmission of start and stop characters along with the barcode value for the Codebar barcode type, an error occurred.</p>
2B0A	<p><b>Priming plate not detected</b></p>
2B0x	<p><b>Dispenser syringe 1 or 2 (respectively) did not home</b></p> <p>x=1–3</p> <p>Generally, this error indicates the syringe was not properly installed. Make sure the syringe's thumbscrews are properly threaded. (Refer to the <b>Installation</b> chapter for instructions.) Restart the reader.</p>
2B04	<p><b>Dispenser syringe 1 or 2 (respectively) failed position verify</b></p> <p>Generally, this error indicates the syringe was not properly installed. Make sure the syringe's thumbscrews are properly threaded. Restart the reader. (Refer to the <b>Installation</b> chapter for instructions.)</p>
37x0/47x0 38x0/48x0 39xy/49xy	<p><b>Noise Test Errors</b> <b>Offset Test Errors</b> <b>Dark Range Errors</b></p> <p>x=0, 1; y=0–6</p> <p>This series of System Test errors may indicate too much light inside the chamber. Make sure the plate carrier door and the front hinged door are properly closed. For models with the dispense module, if the dispense tubes are not connected to the reader, re-install the light shield that shipped with the instrument (or cover the hole with black tape). Restart the reader.</p>

4Fxy	<p><b>Fluorescence signal out of range.</b></p> <p>x=0, 1; y=0–6</p> <p>Too low of a reading indicates a light signal problem. Ensure that Gen5 Fluor/Lum wavelengths table matches the actual filter installed in the filter cube.</p>																													
40xx	<p><b>PMT overload well error at &lt;well #xxx&gt;</b></p> <p>This error typically means that the fluid in a well has oversaturated the PMT (i.e., the well is too bright). Try lowering the sensitivity value in the read step.</p> <p><b>To identify the well:</b> Wells are counted starting at A1, moving left-to-right, row-by-row. The row and column of the well can be extracted from the well number code by applying the following formula (example uses 8 x 12 geometry, 96-well plate):</p> <ol style="list-style-type: none"> <li>1. Convert the ASCII hex string to a decimal equivalent. Ex: “057” indicates 57 hex, yielding a well code of 87 decimal.</li> <li>2. Row = (well code) / (columns in plate), rounded up to a whole number. Ex: 87/12 = 7.25, indicating row 8 (or H).</li> <li>3. Column = (well code) - ((row-1) * (columns in plate)). Ex: 87 - ((8 - 1) * 12) = column 3.</li> </ol> <p>NOTE: If this code is returned during an area scan, it indicates the scan point corresponding to the row/column equivalent in the currently defined scan map, NOT the actual well where the error occurred.</p>																													
4Exy	<p><b>Detector saturated (too much light). Relative Fluorescing Units (RFU) reached (99999).</b></p> <p>X = Fluorescence channel Fluorescence/Luminescence Channel</p> <table border="0" data-bbox="523 1137 1197 1415"> <tr><td>    Mono optics reference channel</td><td>0</td></tr> <tr><td>    Mono optics PMT channel</td><td>1</td></tr> <tr><td>    Top filter optics reference channel</td><td>2</td></tr> <tr><td>    Top filter optics, top PMT channel</td><td>3</td></tr> <tr><td>    Top filter optics, side PMT channel</td><td>4</td></tr> <tr><td>    Bottom filter optics reference channel</td><td>5</td></tr> <tr><td>    Bottom filter optics PMT channel</td><td>6</td></tr> </table> <p>Y = PMT Test Type Code</p> <table border="0" data-bbox="375 1496 1197 1691"> <tr><td>Connection Test</td><td>0</td><td>PMT not connected</td></tr> <tr><td>High Voltage Test</td><td>1</td><td>Failure during test at higher voltage</td></tr> <tr><td>Low Voltage Test</td><td>2</td><td>Failure during test at lower voltage</td></tr> <tr><td>Well Overload Test</td><td>5</td><td>Failure during test at well</td></tr> <tr><td>Background Overload Test</td><td>8</td><td>Failure during background overload test</td></tr> </table> <p>OR Y = filter readset</p> <p>This error can indicate one of several scenarios. It is possibly due to incorrect chemistry, e.g., the fluorescence standards dispensed to the plate exceed expectations.</p> <p>For models with the dispense module, the internal chamber may require cleaning (contact BioTek TAC).</p> <p>If a <b>4E18</b> error is detected during monochromator-based fluorescence, the luminescence probe may be picking up stray light. Try installing a plug in the filter cube. Restart the reader.</p>	Mono optics reference channel	0	Mono optics PMT channel	1	Top filter optics reference channel	2	Top filter optics, top PMT channel	3	Top filter optics, side PMT channel	4	Bottom filter optics reference channel	5	Bottom filter optics PMT channel	6	Connection Test	0	PMT not connected	High Voltage Test	1	Failure during test at higher voltage	Low Voltage Test	2	Failure during test at lower voltage	Well Overload Test	5	Failure during test at well	Background Overload Test	8	Failure during background overload test
Mono optics reference channel	0																													
Mono optics PMT channel	1																													
Top filter optics reference channel	2																													
Top filter optics, top PMT channel	3																													
Top filter optics, side PMT channel	4																													
Bottom filter optics reference channel	5																													
Bottom filter optics PMT channel	6																													
Connection Test	0	PMT not connected																												
High Voltage Test	1	Failure during test at higher voltage																												
Low Voltage Test	2	Failure during test at lower voltage																												
Well Overload Test	5	Failure during test at well																												
Background Overload Test	8	Failure during background overload test																												

	<p><b>Fluorescence wavelength not found in table</b></p> <p>This error indicates that the wavelength specified in the procedure is not detected in the instrument's filter table. In Gen5, verify the Fluorescence filter table has the wavelengths loaded into the reader. Compare the contents of the table with the filters installed in the filter cube (see the Gen5 Help system for more information). Restart the reader.</p>						
50xx 510x	<p><b>Axis failed to home</b></p> <table> <tr> <td>Top filter optics, lower filter/mirror slider</td> <td>03</td> </tr> <tr> <td>Top filter optics, upper filter/mirror slider</td> <td>04</td> </tr> <tr> <td>Bottom filter optics, filter/mirror slider</td> <td>06</td> </tr> </table> <p>Generally, this error indicates the filter cube is not seated properly in the reader. Remove it, ensure each filter or plug is properly positioned and reinstall it securely. Restart the reader.</p>	Top filter optics, lower filter/mirror slider	03	Top filter optics, upper filter/mirror slider	04	Bottom filter optics, filter/mirror slider	06
Top filter optics, lower filter/mirror slider	03						
Top filter optics, upper filter/mirror slider	04						
Bottom filter optics, filter/mirror slider	06						
540x	<p><b>Filter cube failed positional verify</b></p> <table> <tr> <td>Top filter optics, lower filter/mirror slider</td> <td>03</td> </tr> <tr> <td>Top filter optics, upper filter/mirror slider</td> <td>04</td> </tr> <tr> <td>Bottom filter optics, filter/mirror slider</td> <td>06</td> </tr> </table> <p>Generally, this error indicates the filter cube is not seated properly in the reader. Remove it, ensure each filter or plug is properly positioned and reinstall it securely. Restart the reader.</p>	Top filter optics, lower filter/mirror slider	03	Top filter optics, upper filter/mirror slider	04	Bottom filter optics, filter/mirror slider	06
Top filter optics, lower filter/mirror slider	03						
Top filter optics, upper filter/mirror slider	04						
Bottom filter optics, filter/mirror slider	06						

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