


如何导入数据库

Start the *BP-2000 Analysis Software*. 运行 BP-2000分析软件

Click *Data* on the main menu, and then click *Manage Databases*.
点击 main menu, data, 然后 Manage Databases

选择 Import a database



Database Operations


Select Database Operation

Select the database operation you wish to perform:

- Select an existing database
- Create a new database
- Backup (export) a database
- Modify name (or other information) of a database
- Import a database
- Delete a database

Next > Cancel

选择您拷贝的文件夹



Database Operations

Select Import Folder

Select the folder containing the database you wish to import, by typing in the edit box, or using the Browse button.

Source folder:

Browse...

< Back Next > Cancel